



GAMMON TECHNICAL PRODUCTS, INC.  
 P.O. BOX 400 - 2300 HWY 34  
 MANASQUAN, N.J. 08736

PHONE 732-223-4600  
 FAX 732-223-5778  
 EMAIL [gammontech@gammontech.com](mailto:gammontech@gammontech.com)

MATCHED-  
 WEIGHT AND  
 COLOR RATING  
 MEMBRANES

BULLETIN 73  
 (8-03)

PLASTIC MONITORS  
 COLOR RATING MEMBRANES  
 MATCHED-WEIGHT MEMBRANES

APPROVED FOR TESTING JET  
 FUEL PER ASTM METHODS  
 D2276/IP216 AND D5452

Our 37mm, 0.8 µm (micrometer) membranes are produced in exact accordance with ASTM Research Report D1012 under which they have been granted approval for use in methods D2276 as well as IP216. We manufacture our own plastic monitors using Tenite, a product of Eastman Chemical Co. Our monitors are made to fit all recognized brands of holders including the Gammon MiniMonitor Kit<sup>®</sup> and the Millipore<sup>®</sup> Corp. XX64-037-30.

Refill membranes are packaged with a blue paper separator between each one to aid in separation.

Matched-weight refill membranes are packaged with each pair in a glassine envelope. Each pair is matched in weight within 100µg (micrograms) in a closely controlled atmosphere.

Support pads are 34mm in diameter so as to fit into the recess of the plastic monitor beneath the membrane.

Matched -weight monitors have yellow plastic plugs in the inlet port and a red plug in the outlet. Single membrane monitors, used for color rating tests, have blue plugs in the inlet and red plugs in the outlet.

TYPE	MODEL NUMBER	INTERCHANGABLE WITH MILLIPORE P/N	DESCRIPTION ALL MEMBRANES 0.8 MICRON
SINGLE MEMBRANE MONITORS	GTP-1985	MAWP-037-PO	48 PLASTIC MONITORS. EACH WITH A MEMBRANE AND A SUPPORT PAD
MATCHED-WEIGHT MONITORS	GTP-1986	MAWP-037-PM	48 PLASTIC MONITORS. EACH WITH A PAIR OF MEMBRANES AND A SUPPORT PAD
SINGLE MEMBRANE REFILLS	GTP-1983	AAWP-037-PO	100 MEMBRANES AND 100 SUPPORT PADS
MATCHED-WEIGHT REFILLS	GTP-1984	AAWP-037-PM	50 MATCHED-WEIGHT PAIRS OF MEMBRANES AND 50 SUPPORT PADS



MiniMonitor Kit is a registered trademark of Gammon Technical Products. Millipore is a registered trademark of the Millipore Corp.

## MEMBRANE FILING NOTEBOOK AND PAGES

Black 3-ring leatherette notebook holds membrane filing pages. Pages made of heavy manila board, each with pockets to hold 21 mini-envelopes (see below).

NOTEBOOK.....GTP-1400

FILING PAGES.....GTP-1331

packet contains ten (10) sheets - minimum order



## MINI-ENVELOPES

These mini-envelopes are zip-lock polyfilm bags made to insure protection and preservation of filter membranes. Specially prepared white panel for writing data regarding the test such as date, rating, location, etc.

MINI-ENVELOPES...GTP-1267

packet of five hundred (500) envelopes



## TWEEZERS

Our tweezers have blunt points and a slip-lock device. They grip and hold membranes without needing to be held tight in your hand.

TWEEZERS.....GTP-2099



## 47MM MEMBRANES

We have these filter membranes in stock. They are intended for testing jet fuel for particulate content in accordance with ASTM Method D5452. The membranes are white and will fit all brands and models of filter holders that accept 47mm membranes.

### MEMBRANES:

package of 100 membranes, 0.65 $\mu$ m.....TL-3777

package of 100 membranes, 0.8 $\mu$ m.....TL-3777-1